Determine™

HIV-1/2

This package insert must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of assay results cannot be guaranteed if there are deviations from the instructions in this package insert.

NAME AND INTENDED USE

The Abbott Determine™ HIV-1/2 is an *In Vitro*, visually read, qualitative immunoassay for the detection of antibodies to HIV-1 and HIV-2 in human serum, plasma or whole blood. The test is intended as an aid to detect antibodies to HIV-1/HIV-2 from infected individuals.

SUMMARY AND EXPLANATION OF THE TEST

AIDS (Acquired Immunodeficiency Syndrome) is characterized by changes in the population of T-cell lymphocytes. In an infected individual, the virus causes depletion of helper T-cells, which leaves the person susceptible to opportunistic infections and some malignancies. The virus that causes AIDS exists as two related types known as HIV-1 and HIV-2. The presence of the AIDS virus elicits the production of specific antibodies to either HIV-1 or HIV-2.^{1,2,3}

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

Determine HIV-1/2 is an immunochemical test for the qualitative detection of antibodies to HIV-1 and HIV-2.

Sample is added to the sample pad. As the sample migrates through the conjugate pad, it reconstitutes and mixes with the selenium colloid-antigen conjugate. This mixture continues to migrate through the solid phase to the immobilized recombinant antigens and synthetic peptides at the patient window site.

If antibodies to HIV-1 and/or HIV-2 are present in the sample, the antibodies bind to the antigen-selenium colloid and to the antigen at the patient window, forming a red line at the patient window site.

If antibodies to HIV-1 and/or HIV-2 are absent, the antigen-selenium colloid flows past the patient window, and no red line is formed at the patient window site.

To insure assay validity, a procedural control bar is incorporated in the assay device.

CONTENTS

Abbott Determine HIV-1/2 Serum/Plasma Assay (List No. 7D23-13), 100 Tests

- Determine HIV-1/2 Test Card, 10 cards (10 tests/card), HIV-1/2 recombinant antigen and synthetic peptide coated.

Abbott Determine HIV-1/2 Whole Blood Assay (List No. 7D23-33), 100 Tests

- Determine HIV-1/2 Test Card, 10 cards (10 tests/card), HIV-1/2 recombinant antigen and synthetic peptide coated.

- 1 Bottle (2.5 mL) Chase Buffer (List No. 7D22-11) prepared in phosphate buffer. Preservatives: Antimicrobial Agents.

ACCESSORIES (required but not provided)

Serum/Plasma or Whole Blood (venipuncture assay)		Whole Blood (fingerstick assay)	
Pipette	No. 7D22-51	Lancets	No. 7D22-31
Pipette Tips	No. 7D22-61	EDTA Capillary Tubes	No. 7D22-21

WARNINGS AND PRECAUTIONS

For *In Vitro* Diagnostic Use.

CAUTION:

Appropriate biosafety practices^{4,5} should be used for materials that contain or are suspected of containing infectious agents. These precautions include, but are not limited to the following:

- Wear gloves when handling specimens or reagents.
- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect all spills of specimens or reagents using a tuberculocidal disinfectant, such as 0.5% sodium hypochlorite, or other suitable disinfectant.^{6,7}
- Decontaminate and dispose of all specimens, reagents, and other potentially contaminated materials in accordance with local regulations.^{8,9}
- Kit components are stable until expiration date when handled and stored as directed. Do not use kit components beyond expiration date.

STORAGE

The Abbott Determine HIV-1/2 Test Cards and Chase Buffer must be stored at 2-30°C until expiration date.

SPECIMEN COLLECTION

Serum, Plasma, and Whole Blood Collection by Venipuncture

Human serum, plasma, and whole blood collected by venipuncture should be collected aseptically in such a way as to avoid hemolysis.

NOTE: For whole blood and plasma specimens, EDTA collection tubes must be used.

Whole Blood Collection by Fingerstick¹⁰

Before collecting a fingerstick specimen, place an EDTA capillary tube on a clean dry surface.

- Choose the fingertip of the middle, ring, or index finger (whichever is the least calloused) for adults and children older than one year. Warm the hand as needed with a warm, moist towel or warm water to increase blood flow.

- Clean fingertip with alcohol; allow to air dry.

- Position the hand palm-side up. Place the lancet off-center on the fingertip. Firmly press the lancet against the finger and puncture the skin. Dispose of the lancet in an appropriate biohazard sharps container.

- Wipe away the first drop of blood with a sterile gauze pad.

- Hold the finger lower than the elbow and apply gentle, intermittent pressure to the base of the punctured finger several times. Touch the tip of the EDTA Capillary Tube to the drop of blood*. Avoid air bubbles.

*If EDTA Capillary Tubes (No. 7D22-21) will be used, fill the tube with blood between the 2 marked lines.

SPECIMEN STORAGE

- Serum and plasma specimens should be stored at 2-8°C if the test is to be run within 7 days of collection. If testing is delayed more than 7 days, the specimen should be frozen (-20°C or colder).
- Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 7 days of collection. Do not freeze whole blood specimens.
- Whole blood collected by fingerstick should be tested immediately.

TEST PROCEDURE

The desired number of test units from the 10-test card can be removed by bending and tearing at the perforation.

NOTE: Removal of the test units should start from the right side of the test card to preserve the lot number which appears on the left side of the test card.

- Remove the protective foil cover from each test.
- For serum or plasma samples:
 - Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
 - Wait a minimum of 15 minutes (up to 60 minutes) and read result.
- For whole blood (venipuncture) samples:
 - Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
 - Wait one minute, then apply one drop of Chase Buffer to the sample pad.
 - Wait a minimum of 15 minutes (up to 60 minutes) and read result.
- For whole blood (fingerstick) samples:
 - Apply 50 µL of sample (by EDTA capillary tube) to the sample pad (marked by the arrow symbol).
 - Wait until blood is absorbed into the sample pad, then apply one drop of Chase Buffer to the sample pad.
 - Wait a minimum of 15 minutes (up to 60 minutes) and read result.

QUALITY CONTROL

To insure assay validity, a procedural control is incorporated in the device and is labeled "Control". If the control bar does not turn red by assay completion, the test result is invalid and the sample should be retested.

INTERPRETATION OF RESULTS

POSITIVE (Two Bars)

Red bars appear in both the control window (labeled "Control") and the patient window (labeled "Patient") of the strip. Any visible red color in the patient window should be interpreted as positive.


NEGATIVE (One Bar)

One red bar appears in the control window of the strip (labeled "Control"), and no red bar appears in the patient window of the strip (labeled "Patient").

INVALID (No Bar)

If there is no red bar in the control window of the strip, and even if a red bar appears in the patient window of the strip, the result is invalid and should be repeated. If the problem persists, contact your local Abbott Customer Service and Support Center.

NOTES:

- The test result is positive even if the patient bar appears lighter or darker than the control bar.
- If an invalid test result occurs repeatedly, or for technical assistance, contact your local Abbott Customer Service and Support Center.

LIMITATIONS OF THE PROCEDURE

- The Abbott Determine HIV-1/2 test is designed to detect antibodies to HIV-1 and HIV-2 in human serum, plasma, and whole blood. Other body fluids or pooled specimens may not give accurate results.
- The intensity of the patient bar does not necessarily correlate to the titer of antibody in the specimen.
- No test provides absolute assurance that a sample does not contain low levels of antibodies to HIV-1 and HIV-2 such as those present at a very early stage of infection. Therefore a negative result at any time does not preclude the possibility of exposure to or infection with HIV-1 or HIV-2 viruses.
- Positive specimens should be retested using another method and the results should be evaluated in light of the overall clinical evaluation before a diagnosis is made.
- Whole blood or plasma specimens containing anticoagulants other than EDTA may give incorrect results.

PERFORMANCE CHARACTERISTICS

SPECIFICITY

A total of 1,594 serum and plasma specimens from Asia, West Africa, and North America were tested by Abbott Determine HIV-1/2 and a commercially available test (Table I).

Table I
Specificity of Abbott Determine HIV-1/2

Population	Number of Specimens Tested	Negative by Abbott Determine HIV-1/2	Negative by a Commercially Available Test
Seronegative Serum Plasma	908 403	907/908 (99.89%) 403/403 (100.00%)	908/908 (100.00%) 403/403 (100.00%)
Pregnant Females	58*	57/57 (100.00%)	57/57 (100.00%)
West Africans	49	48/49 (97.96%)	48/49 (97.96%)
Disease States Other than HIV and Potentially Interfering Substances	176*	173/175 (98.86%)	174/175 (99.45%)
Total	1,594	1,588/1,592 (99.75%)	1,590/1,592 (99.87%)

*One specimen from a pregnant female and an HCV positive patient were positive by both Abbott Determine and the commercially available test. Both specimens confirmed positive by HIV-1 Western Blot.

A total of 368 seronegative whole blood specimens from Thailand were tested with paired serum and plasma by Abbott Determine HIV-1/2. Thirty-nine of the whole blood specimens were collected by both venipuncture and fingerstick (Table II).

Table II
A Comparison of Abbott Determine HIV-1/2 Specificity in Seronegative Whole Blood and Paired Serum and Plasma Specimens

Specimen Type	Number of Specimens Tested	Negative by Abbott Determine HIV-1/2
Serum	368	368/368 (100.00%)
Plasma	368	368/368 (100.00%)
Whole Blood (venipuncture)	368	368/368 (100.00%)
Whole Blood (fingerstick)	39	39/39 (100.00%)

SENSITIVITY

A total of 869 HIV-1 and HIV-2 antibody positive serum and plasma specimens from Asia, Africa, North and South America were tested by Abbott Determine HIV-1/2 and a commercially available test (Table III).

Table III
Sensitivity of Abbott Determine HIV-1/2

Population	Number of Specimens Tested	Positive by Abbott Determine HIV-1/2	Positive by a Commercially Available Test
HIV-1 Positive	521	521/521 (100.00%)	521/521 (100.00%)
HIV-2 Positive	114	114/114 (100.00%)	114/114 (100.00%)
HIV-1 Subtypes A-G	222	222/222 (100.00%)	Not Tested Not Tested
HIV-1 Group O	12	12/12 (100.00%)	Not Tested Not Tested
Total	869	869/869 (100.00%)	635/635 (100.00%)

A total of 102 seropositive whole blood specimens from Thailand were tested with paired serum and plasma by Abbott Determine HIV-1/2. Thirty-two of the whole blood specimens were collected by both venipuncture and fingerstick (Table IV).

Table IV
A Comparison of Abbott Determine HIV-1/2 Sensitivity in Seropositive Whole Blood and Paired Serum and Plasma Specimens

Specimen Type	Number of Specimens Tested	Positive by Abbott Determine HIV-1/2
Serum	102	102/102 (100.00%)
Plasma	102	102/102 (100.00%)
Whole Blood (venipuncture)	102	102/102 (100.00%)
Whole Blood (fingerstick)	32	32/32 (100.00%)

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BIBLIOGRAPHY / LITERATUR / BIBLIOGRAPHIE / BIBLIOGRAFÍA / BIBLIOGRAFÍA

- Piot P, Plummer FA, Mhalu FS, Lambary JL, Chin J, Mann JM. AIDS: An International Perspective. *Science*. 1988;239:573-579.
 - Weniger BG, Takebe Y, Ou CY, et al. The Molecular Epidemiology of HIV in Asia. *AIDS*. 1994;8(S2):S13-S28.
 - Günler LG, Hauser PH, Eberle J, et al. A New Subtype of Human Immunodeficiency Virus Type 1 (MVP-5180) from Cameroon. *Journal of Virology*. 1994;68(3):1581-1585.
 - World Health Organization. *Laboratory Biosafety Manual*. Geneva: World Health Organization,1993.
 - National Committee for Clinical Laboratory Standards. *Protection of Laboratory Workers from Infectious Disease Transmitted by Blood, Body Fluids, and Tissue*. Tentative Guideline. NCCLS Document M29-T2. Villanova, PA: NCCLS, 1991:1-43.
 - CDC. Recommendations for Prevention of HIV Transmission in Health-Care Settings, *MMWR* 1987;36(25):3S-18S.
 - Schulster LM, Hollinger FB, Dreesman GR, et al. Immunological and Biophysical Alteration of Hepatitis B Virus Antigens by Sodium Hypochlorite Disinfection. *Appl Environ Microbiol* 1981;42(S):762-7.
 - National Committee for Clinical Laboratory Standards. *Clinical Laboratory Waste Management*: Approved Guideline. NCCLS Document GP5-A. Villanova, PA: NCCLS, 1993;13(22):1-18, 29-42.
 - US Environmental Protection Agency. *EPA Guide for Infectious Waste Management*. Publication No. EPA/530-SW-86-014. Washington, DC: US Environmental Protection Agency, 1986;1-1 p-5, R1-R3, A1-A24.
 - NCCLS Guidelines, No. H4A3, Volume 11, No. 11, Approved Standard, July 1991.

Deutschland	France	México
ABBOTT GmbH Diagnostika Max-Planck-Ring 2 D-65205 Wiesbaden-Delkenheim Tel. (061 22) 580 Fax (06122) 58 12 44	Laboratoires ABBOTT Division Diagnostic Tissue et la Couture SILIC 203 F-94518 RUNGIS Cedex Tel. (01) 45 60 25 00 Fax (01) 45 60 04 98	ABBOTT Laboratorios de México S.A. de C.V. Apartado Postal No. 44-983 México 03100, D.F. Tel. (525) 726-4600 Fax (525) 726-4644
Australia	Great Britain	Puerto Rico
ABBOTT Australia Pty Ltd. Diagnostics Division P.O. Box 394 North Ryde N S W 2113 Tel. (02) 888-1111 Fax (02) 887-3948	ABBOTT DIAGNOSTICS Abbott House Norden Road Maidenhead Berkshire SL6 4XL Tel. (06 28) 78 40 41 Fax (06 28) 64 62 05	ABBOTT Laboratories Inc. Diagnostics Division P.O. Box 4706 Puerto Rico 00628 Tel. (787) 762-3368 Fax (787) 750-4428
Belgium / Luxembourg	Hellas	Schweiz
ABBOTT S.A./N.V. Diagnostics Division Rue du Bosquet 2 B-1348 Ottignies Louvain-la-Neuve Tel. (010) 47 53 11 Fax (010) / 47 53 34	ABBOTT Laboratories (Hellas) S.A. Diagnostic Division 194 Syngrou Avenue GR-176 71 Kallithea Athens Tel. (01) 95 19-019 Fax (01) 95 92-790	ABBOTT A.G., Diagnostics Division Gewerbestrasse 5 CH-6330 Cham Tel. (041) / 74 94 422 Fax (041) / 74 11 802
Canada	Hong Kong	Singapore
ABBOTT Laboratories, Ltd. Diagnostics Division 777 North Service Creek Drive Mississauga, Ontario LGN 3R3 Tel. (416) 858-2450 Fax (416) 858-2462	ABBOTT Laboratories Ltd. Greater China & Hong Kong B 605, Seaview Estate 2-8 Watson Road North Point, Hong Kong Tel. (852) 566-8711 Fax (852) 806-3922	ABBOTT Laboratories Private Ltd. P.O. Box 1016 Singapore 9020 Tel. 278-8343 Fax 270-8873
Danmark	Ireland	South Africa
ABBOTT Laboratories A/S Diagnostics Division Bygstrupben 16 Troerod DK-2950 Vedbaek Tel. (45) 45 67 02 00 Fax (45) 45 67 02 02	ABBOTT Laboratories Ireland Limited 70, Broomhill Road Dublin 24 Tel. (016) 51 73 88 Fax (016) 51 77 65	ABBOTT Diagnostics Division P.O. Box 1616 Johannesburg 2000, South Africa Tel. (011) 494-7111 Fax (011) 494-7101
España	Italia	Norway
ABBOTT CIENTIFICA, S.A. División Diagnósticos Costa Brava, 13 28034 Madrid Tel. (91) 337 34 00 Fax (91) 734 96 64	ABBOTT S.p.A. Divisione Diagnosi Via Mar della Cina 262 I-00144 Roma (RM) Tel. (06) 529-911 Fax (06) 529-914 36	ABBOTT NORGE Diagnostic Division Postboks 822 Nesoytjern 4 N-1361 Billingstad Tel. (47) 68 84 55 00 Fax (47) 68 84 84 43
Finnland	Portugal	Sverige
ABBOTT Diagnostics Vapaalahti 2A SF-01650 Vantaa Tel. (0) 653-4788 Fax (0) 653-1815	Dainabot Co., Ltd. Roppongi First Building 9-9, Roppongi 1-chome Minato-ku, Tokyo 106 Tel. (03) 35 89-9441 Fax (03) 35 89-9591	ABBOTT Scandinavia AB Diagnostics Division P.O. Box 1074 S-16421 Kista-Stockholm Tel. (08) 703-6700 Fax (08) 751-3060
Japan	Philippines	USA
	ABBOTT Laboratories Diagnósticos 2700 Amadora Tel. (01) 471-2725 Fax (01) 471-2725	ABBOTT Laboratories Diagnosis Division 100 Abbott Park Road Abbott Park, IL 60064-3500 Tel. (847) 937-6100
Portugal	Philippines	ABBOTT Laboratories
	ABBOTT Laboratories Diagnósticos 2700 Amadora Tel. (01) 471-2725 Fax (01) 471-2725	Diagnostics Division 820 Mission Street South Pasadena, CA 91030 Tel. (818) 440-0700
Puerto Rico	Philippines	ABBOTT Laboratories
	ABBOTT Laboratories Diagnósticos 2700 Amadora Tel. (01) 471-2725 Fax (01) 471-2725	Diagnostics Division 1921 Hurst Drive Irving, TX 75038 Tel. (214) 518-6000