



# *Post-marketing Activities*

## *Industry Perspective*

Luann Ochs, MS  
CLIAC, Feb 11, 2004





# The Regulations



## *Quality System Regulation, 21 CFR Part 820*

- 820.198 Complaint Files

Each manufacturer shall establish procedures for receiving, reviewing, and evaluating complaints

- Timely processing
- Documented
- Determine whether the complaint represents an event that is required to be reported to FDA

## *Quality System Regulation, 21 CFR Part 820*

- 820.100 Corrective and Preventive Action (CAPA)

Each manufacturer shall establish and maintain procedures for implementing corrective and preventive action

- Analyze processes, quality audit reports, quality records,

- service records, complaints, returned product

- Identify existing and potential causes of nonconforming

- product or quality problems

- Identify actions needed to correct and prevent recurrence

- Verify or validate corrective or preventive action

# The Regulations, continued



## *Medical Device Reporting (MDR), 21 CFR 803*

Device user facilities, importers, and manufacturers must report to FDA when a device has, or may have caused or contributed to a death or serious injury

- Adverse event files
- Submit summary reports to FDA
- Time requirements: 30 days, or 5 days if remedial action is needed to prevent further risk of harm to public health

## *Companies utilize customer calls for*

- Fulfilling Regulatory obligations, as well as for
- Continuous product improvement
- Promote customer satisfaction



# Customer Calls the Hotline

## *Allegation of death or serious injury?*

### YES

- Get all details of incident
- Follow procedures for MDR reporting to FDA
- Analyze retentions or returns, looking for product quality issues
- CAPA?

### NO


- Phone solve
- Troubleshoot
  - May send product
- Replace
  - May request product return for analysis
- Trend analysis
- CAPA?



Risk

# Trend Analysis

- Vast majority of CAPAs are result of trend analysis
- Complaint data compiled over time
- Break down data, e.g.:
  - by product
  - by specific complaint
  - by resolution of complaint
  - by identified failure






# Total Calls & Resolution




**OTC / CLIA-waived Glucose Meter**  
**Quarterly Call Summary**



# Returned Product Analysis




**OTC / CLIA-waived Glucose Meter**  
**Analysis of Returned Product by Complaint**



# Return Analysis, Normalized




**OTC / CLIA-waived Glucose Meter**  
Complaints normalized for production volume



# Failure Analysis



**OTC / CLIA-waived Meter  
12 Month Component Failure Analysis**



## 1. *Identify the hazards*

- FMEA, FTA  
Example FMEA template
- NCCLS EP 18  
EP-18 Examples



*Source: FDA Guidance "Device Use Safety, Incorporating Human Factors in Risk Management, Aug 1999.*

# Example FMEA template



Possible Hazard	Effect	Potential Causes	Minimum Requirement for safety	Controls	S	O	D	RPN

## 1. *Identify the hazards*

- FMEA, FTA  
Example FMEA template
- NCCLS EP 18  
EP-18 Examples



# EP-18 Example 1



**Appendix A. Example of a “System-Specific Sources of Error” Matrix**

Potential Sources of Error	Applicable Y/N?	Nature of Impact	Training/Laboratory Procedure Requirements	Applicable Quality Monitoring	Frequency of Monitoring
<b>1 Specimen Collection</b>					
1.1 Contamination					
1.1.1 Alcohol					
1.1.2 Other Cleansing Agent					
1.1.3 Anticoagulants in Lines					
1.1.4 Intravenous Fluids					
1.1.5 Admixture with Other Fluids/Materials					
1.2 Inadequate Sample					
1.2.1 Poor Circulation at Sample Site					
1.2.2 Poor Vascular Access					
1.2.3 Not Enough Collected					
1.2.4 Poor Technique					
1.2.5 Too Much Collected					
1.3 Hemolysis					
1.4 Incorrect Patient Drawn					
1.5 Inappropriate Sample					
1.5.1 Arterial vs. Venous vs. Capillary					
1.5.2 Whole Blood vs. Plasma					
1.5.3 Sample in Wrong Container or Syringe/Wrong Additives					
1.5.4 Fasting vs. Nonfasting					
1.5.5 Clotted Sample					
1.5.6 Inappropriate Time of Collection					
1.6 Patient Condition Inappropriate for Testing Method					



# EP-18 Example 2



## Appendix A. (Continued)

Potential Sources of Error	Applicable Y/N?	Nature of Impact	Training/Laboratory Procedure Requirements	Applicable Quality Monitoring	Frequency of Monitoring
4.5 Calculation Required					
<b>5 Preliminary Review</b>					
5.1 Improper Interpretation of Control Results					
5.2 Outlier/Nonsense Result Not Recognized					
5.3 Result Outside of Linear Range Not Recognized					
5.4 Alert Value Not Recognized					
5.5 Need for a Confirmatory Sample Not Recognized					
5.6 Effect of Preanalytical Variables Not Recognized					
5.7 Instrument Malfunction Not Recognized					
5.8 Interference Not Recognized					
<b>6 Integration/Report into Chart</b>					
6.1 No Result Recorded					
6.2 Result Recorded in Incorrect Patient Chart					
6.3 Incorrect Information Recorded					
6.3.1 Data					
6.3.2 Time					
6.3.3 Result					
6.4 Information Unreadable					

# Risk Management Principles, continued

## 2. *Develop strategies and controls to eliminate, reduce likelihood of, or mitigate the consequences*

- Remove hazard causes through design
- Make design error tolerant
- Alert users
- Develop written procedures and training for safe operation (labeling and training materials)



# Risk Management Principles, continued

## 3. *Verify that controls are effective,* *e.g.:*

- Validation
- Stress testing
- Human factors testing
- User studies



# Summary

