

MATERIAL SAFETY DATA SHEET

1. COMPANY IDENTIFICATION

BIOKIT S.A.
08186 Llíssà d'Amunt (Barcelona)
Spain
Telephone: +34-93-860 9000
Fax : +34-93-860 9009

2. PRODUCT NAME

biorapid HIV 1+2

Code # 3000-1069

3. HAZARDOUS COMPONENTS

Not apply.

4. INTENDED USE

Rapid qualitative two step immunoassay test for the detection of antibodies to HIV-1 and HIV-2 in human serum, plasma or whole blood.

5. PHYSICAL AND CHEMICAL PROPERTIES

Not apply.

6. HAZARDS TO HEALTH

When used in accordance with the Principles of Good Laboratory Practice, Good Standards of Occupational Hygiene and the instructions stated in the product insert, these products are not considered to present a hazard to health.

7. PRECAUTIONS

The storage must be between 4 and 30°C.

For IN VITRO diagnostic use.

8. EMERGENCY ACTION

Not apply.

9. FIRE HAZARDS

Extinguishing media:

Carbon dioxide, dry chemical powder or polymer foam.

Special firefighting procedures:

wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.

10. DISPOSAL REQUIREMENTS

Remains of samples, reagents and controls should be collected in a recipient for this purpose and autoclaved 1 hour at 121°C.

Observe all Federal, State and local laws.

11. TRANSPORT PRECAUTIONS

Special requirements: NONE.

THE ABOVE INFORMATION IS BELIEVED TO BE CORRECT, BUT DOES NOT PURPOSE TO BE ALL INCLUSIVE AND SHOULD BE USED ONLY AS A GUIDE. biokit SHALL NOT BE HELD FOR ANY DAMAGE RESULTING FROM HANDLING OR USE OF THE PRODUCT.

Quality Assurance Department.