
# Antimicrobial Susceptibility Testing (AST): Current Issues & Implications for Public Health

Janet Fick Hindler, MCLS MT(ASCP)
jhindler@ucla.edu
UCLA Medical Center
Los Angeles, CA



# At the conclusion of this talk, you will be able to.....

- Discuss current issues in antimicrobial susceptibility testing and antimicrobial resistance.
- List current sources for information on antimicrobial susceptibility testing.
- Describe the effects of these issues on public health.


#### **Current Issues in AST**

- Identify methods to detect and report "emerging resistance"
  - Reliable
  - Practical for clinical laboratories
- Communicate AST and reporting issues to clinical labs
- Assess performance of clinical labs in detecting "emerging resistance"

# Communicate AST & Reporting Issues to Clinical Labs

- Primary resources for continuing education
  - Workshops / seminars
  - Internet
  - Publications
- CDC / NLTN efforts
  - MASTER website (www.phppo.cdc.gov/dls/master/default.asp)
  - AST CD ROM (available late 2002)
  - NLTN workshops, teleconferences, other

### Detect / Report Emerging Resistance



Clinical Lab

Public Health Lab (Define role)

# Public Health Lab as A Resource for Clinical Labs in AST & Reporting

- Performs AST (confirms unusual results)
- Answers questions
  - Identifies other resources, when needed
- Provides news "bulletins"
- Provides CE

Much Variability Among PH Labs!

### Clinical Lab Need

Knowledge of "resources" for AST & reporting in their community ....particularly when a "real problem" arises.

### Assess Performance of Clinical Labs in Detecting Emerging Resistance

- Inspections
- Proficiency surveys
- Post CE monitoring
- (Competency assessment)



# Contemporary Resistance Concerns

- Staphylococcus aureus
  - -MRSA new tests
  - -VISA
  - -VRSA
- Streptococcus pneumoniae
  - -Reporting using new breakpoints

# Contemporary Resistance Concerns (con't)

- Streptococcus spp.
  - Beta Streptococcus Group B (in pregnancy) and clindamycin / erythromycin
- Enterobacteriaceae
  - Extended-spectrum beta-lactamases (ESBLs)

#### NCCLS Standards - 2002

- M2-A7 Disk Diffusion
- ♦ M7-A5 MIC
- M100-S12 Tables\*

\*Updated annually www.nccls.org (610) 688-0100



# Staphylococcus spp.

Organism	1st Choice Drugs	Alternative Drugs
MSS	P'ase stable penicillin (e.g., oxacillin)	a cephalosporin, vancomycin, ß-lac/ß-lac inhibitor combo, imipenem or meropenem, clindamycin, a fluoroquinolone
MRS	vancomycin +/- gentamicin +/- rifampin	linezolid, quin-dalfo, a fluoroquinolone, a tetracycline, trim-sulfa

# MRSA (ORSA)

- Genetic determinant mecA
- ◆ MecA codes for PBP2a




#### **MRSA Test Methods**

- NCCLS disk diffusion or MIC
- Commercial method
- Agar screen
- mecA assay
- ♦ PBP2a assay

### **MRSA**

"Isolates of staphylococci that are shown to carry the *mecA* gene, or that produce PBP2a, the gene product, should be reported as oxacillin resistant."





PBP2a



# Can mecA or PBP2a Assays be Used Alone?

- MRSA?
- MSSA?
- Are these tests consistently reliable "on the bench"?
- Are results for other drugs routinely needed?

# Staphylococcus - ß-Lactams

Ox	Comments
S	S to penicillins, cephems,
	carbapenems
S	R to ß-lactamase labile pens; S to ß-lactamase stable pens;
	S to ß-lac / ß-lac inhibitor
	combos, cephems, carbapenems
R	R to all ß-lactams
	S

# Staphylococcus aureus

clindamycin S

erythromycin S

oxacillin S

penicillin R

vancomycin S

"Cefazolin and other beta-lactams (except amoxicillin, ampicillin, and penicillins) are active against oxacillin-S and penicillin-R staphylococci."

# VISA\* (GISA\*\*)

- Mostly MRSA
- ♦ Vanco MIC=8 µg/ml
- Japan, USA (MI, NY, NJ, CA), Europe
- Patients previously Rx with vanco
  - \* vancomycin-intermediate *S. aureus*\*\*glycopeptide-intermediate *S. aureus*(vancomycin is a glycopeptide)

#### **Detection of VISA**

MicroScan ON

Etest

Sensititre

Vitek

◆ BHI-Van (6 µg/ml)

MicroScan rapid

Disk diffusion

MIC (µg/ml)

8

6-8

4,8

4

growth


≤**2**, ≥ **16** 

inadequate

Tenover, et al. 1998. JCM. 36:1020


# Case Study VISA - Pt. JB (UCLA)

- ♦27 y.o. referral patient
- liver transplant candidate
- Bile drainage (liver abscesses)
- 2 strains VISA



# Case Study - VISA Pt. JB

Strain	mecA	MIC ()	ug/ml) Van
1	_	0.5	8
2	+	>16	8



### Confirmation of VISA

- If VISA suspected:
  - 1. Repeat ID and susceptibility tests
  - 2. Contact your institution's infection control department
  - 3. Contact CDC at SEARCH@cdc.gov
  - 4. Contact your local health department
  - 5. Save isolate

#### VRSA – 1<sup>st</sup> Case

- Michigan
- 40 y.o., diabetic on dialysis
- Previous catheter site with MRSA
- Previous vancomycin Rx (6.5 wks)

# VRSA - 1st Case (con't)

- ◆VRSA with mecA and vanA genes
- Patient also had VRE (vanA)
- ◆Vancomycin MIC =1024 µg/ml
- S to chloramphenicol, linezolid, minocycline, quin-dalfo, tetracycline, TMP-SMZ

MMWR 51:565-567, July 5, 2002

#### VRSA - 2<sup>nd</sup> Case

- Pennsylvania
- Chronic foot ulcer
- VRSA with mecA and vanA genes
- Vancomycin MIC = 32 μg/ml
- ◆ DD zone = 12 mm
- Van (6 μg/ml) BHI screen = growth
- S to chloramphenicol, linezolid, minocycline, quin-dalfo, rifampin, TMP-SMZ

MMWR 51:902. Oct. 11, 2002

If pen  $\leq$  0.06 µg/ml (S)...

**1st Choice Drugs** 

**Alternative Drugs** 

penicillin G or V; amoxicillin a cephalosporin; erythromycin; clarithromycin; azithromycin; levofloxacin, gatifloxacin, or moxifloxacin; meropenem; imipenem; trimsulfa; clindamycin; a tetracycline

If pen 0.12 - 1 µg/ml (I)...

**1st Choice Drugs** 

penicillin G IV or ceftriaxone or cefotaxime

#### **Alternative Drugs**

levofloxacin, gatifloxacin, or moxifloxacin; vancomycin clindamycin

If pen ≥ 2 μg/ml (R)...

**1st Choice Drugs** 

**Alternative Drugs** 

**Meningitis:** 

vancomycin +

cefotaxime or ceftriaxone

+/- rifampin

Other:

as above or levofloxacin, gatifloxacin, or moxifloxacin

meropenem; imipenem

quinupristin/dalfopristin linezolid

# Streptococcus pneumoniae Breakpoints

#### Ceftriaxone and Cefotaxime

- Originally based on treatment of meningitis
- NCCLS saw need for breakpoints for treatment of respiratory infections
- New non-meningitis breakpoints

#### Penicillin

- Originally based on treatment of meningitis
- No change in breakpoints

MIC (μg/ml) S Int R

Ceftriaxone or cefotaxime (meningitis)

**≤ 0.5** 

1

**≥ 2** 

Ceftriaxone or cefotaxime (non-meningitis)

**≤ 1** 

2

**≥ 4** 

Penicillin

**≤ 0.06** 

0.12 - 1

≥ 2

# Streptococcus pneumoniae Penicillin

"High doses of IV penicillin (e.g. at least 2 million units every 4 hours in adults with normal renal function ) or similarly ampicillin (e.g. 2 g at every 6 hours) are effective in treating pneumococcal pneumonia due to strains in the intermediate category."

### S. pneumoniae (CSF)

MIC (µg/ml)

ceftriaxone

(meningitis) ≤0.25 S

penicillin ≤0.03 S

vancomycin 0.5 S

Patients with meningitis require therapy with maximum doses of ceftriaxone.

### S. pneumoniae (blood)

MIC (µg/ml)
1 I
1 S
R
0.5 S
1.0 I
0.5 S

Pts. with meningitis require therapy with max doses of ceftriax; High dose IV pens (e.g. at least 2 mil U every 4 h in adults with normal renal function) or amp (e.g. 2 g at every 6 h) are effective in treating pneumococcal pneumonia due to strains in the penicillin "int" category.

## Beta Streptococcus spp.

	1st Choice Drugs	Alternative Drugs
Beta Group A, C, G	penicillin G or V	clindamycin; erythromycin; a cephalosporin; vancomycin; clarithromycin; azithromycin
Group B	penicillin G or ampicillin	a cephalosporin; vancomycin; erythromycin

The Medical Letter; 2001; 43:69.

## **Throat Culture**

Many Group A Streptococcus

"Group A Streptococcus remains universally susceptible to penicillin"

# Beta Streptococcus spp. Erythromycin/Clindamycin

- Erythromycin and/or clindamycin are alternatives in penicillin-allergic patients
- Group A 5-10% erythromycin-R
- Group B up to 25% erythromycin-R and 15% clindamycin-R (prophylaxis in pregnant women issue)

# Extended-Spectrum ß-lactamases (ESBLs)

- Result from mutation of common ßlactamase genes (e.g., bla<sub>TEM-1</sub>, bla<sub>SHV-1</sub>)
- Inactivate extended-spectrum ß-lactam agents
- Genes located on plasmids (often with other R genes)
- Approx. 100 TEM and 30 SHV types

## Significance of ESBLs

- Cause nosocomial infections
- Serious, infections. Treatment failures (esp. bacteremia) when patient treated with 3rd gen. cephalosporins.

## ESBLs — Testing

- Organisms:
   E. coli, Klebsiella spp.
- Screen test:
   Decreased susceptibility to extended-spectrum
   ß-lactams
- Phenotypic confirmatory test:
   ß-lactam activity restored by ß-lactamase inhibitor (e.g. clavulanic acid)
- Other clues:
  "S" to cephamycins; multi-R

## Klebsiella pneumoniae ?ESBL

PRELIM:		MIC (µg/ml)
C	efoxitin	1 S
Ci	profloxacin	0.5 S
in	nipenem	≤ <b>0.25 S</b>
pi	p-tazo	8 S
aı	m, cfaz, gm, T-S	R
	***hold cofonime	ofotovimo if "C"

\*\*\*hold cefepime, cefotaxime, if "S"

This *K. pneumoniae* is suspicious for extended-spectrum betalactamase (ESBL) production; confirmatory tests pending.

# **ESBL**Phenotypic Confirmatory Test

#### Test:

- cefotaxime
- cefotaxime/clavulanic acid
- ceftazidime
- ceftazidime/clavulanic acid

#### Results:

clavulanic acid restores activity of cefotaxime or ceftazidime or both

# **ESBL Confirmatory Test**



## **ESBL** Reporting

- Confirmed ESBLs report "R" for:
  - -cephalosporins (not cephamycins)
  - -penicillins (not ß-lac inhibitor combos)
  - -aztreonam

### Klebsiella pneumoniae ESBL

<b>FINAL:</b>		MIC (µg/ml)
	cefepime	R
	cefoxitin	1 S
	cefotaxime	<b>♥</b> R
	ciprofloxacin	0.5 S
	imipenem	0.25 S
	pip-tazo	8 S
	am. cefaz. gent. T-S	R

Confirmatory tests for this *K. pneumoniae* indicate unusual resistance [extended-spectrum beta-lactamase (ESBL) production]; ID consult suggested.

### Issues - ESBLs

- Reporting ESBLs from urine isolates (e.g. acute cystitis)
- Changing results for other β-lactams (e.g. piptazobactam)
- Detection of.....
  - all types of ESBLs (≅100 TEM and 30 SHV)
  - ESBLs present with other "R" mechanisms
  - ESBLs beyond E. coli and Klebsiella spp.
  - other  $\beta$ -lactamases (e.g. plasmid-mediated ampC)

## **NCCLS M39-A Guideline**

"Analysis and Presentation of Cumulative Antimicrobial Susceptibility Test Data"

Focus - clinical laboratories
Guide MDs on empiric therapy



# % Susceptible 2001 – Exmp.

	n	Am	Cf	Ctx	Cip	Gm	Pp	T-S
E. coli	729	66	81	98	97	97	66	76
E. cloacae	144	-	-	66	95	88	65	88
P. aerug	221	-	-	10	79	81	76	-

# **CAP Checklist MIC.21950 (2001)**

- Does the procedure manual address unusual or inconsistent antimicrobial susceptibility testing results?
  - Testing QC strains doesn't guarantee accurate patient results
  - Unusual or inconsistent results should be investigated

# Verification Rules - Exmps. Biologically implausible or infrequent

- Repeat test
  - -amik-R + gent-R + tob-R Enterobacteriaceae
  - -amp-R E. faecalis
  - imipenem-R Enterobacteriaceae
- Reexamine or repeat test
  - -amp-S Klebsiella spp.

## Enterobacter aerogenes

ampicillin R

cefazolin R

ceftizoxime S

gentamicin S

imipenem R\*

trimeth-sulfa S

\* "R" likely due to drug deterioration

## Conclusions

#### Clinical laboratory detection of emerging resistance

- Labs must have access to current testing / reporting guidelines
- Lab personnel must be competent in testing
- Resources for "problems" in emerging resistance must be available to labs

#### Staphylococcus aureus

- mecA or PBP2a tests can be used to identify MRSA
- VISA are uncommon and may be difficult to detect
- VRSA are rare and can be detected readily; however the resistance may be unstable

# Conclusions (con't)

#### Streptococcus pneumoniae

- The new NCCLS breakpoints for ceftriaxone and cefotaxime should be reported
- MICs should be performed on CSF isolates ASAP
- A comment should be added to explain the meaning of a penicillin "Int" result on non-CSF isolates

#### Streptococcus spp.

 Group B streptococci from penicillin-allergic pregnant women should be tested with clindamycin and/or erythromycin

## Conclusions (con't)

#### Enterobacteriaceae

- E. coli and Klebsiella spp. should be tested for ESBL production
- Urine isolates associated with acute uncomplicated cystitis may not require ESBL testing
- There are currently no standard NCCLS methods for detecting ESBLs in isolates other than *E. coli* and *Klebsiella* spp.

#### Local resistance patterns

 When possible, NCCLS M39-A should be followed for preparing a cumulative antibiogram report

# Conclusions (con't)

#### Verifying results on patient isolates

 Each laboratory should develop a written protocol for verification of unusual or atypical patient results prior to reporting them

#### Resources

 a listing of antimicrobial susceptibility testing resources can be found on MASTER website

### **Additional References**

Livermore, D. M., T. G. Winstanley, and K. P. Shannon. 2001. Interpretative reading: recognizing the unusual and inferring resistance mechanisms from resistance phenotypes. J. Antimicrob. Chemother. 48:87-102.

#### Websites:

http://www.cdc.gov/drugresistance/factsheets/index.htm CDC antimicrobial susceptibility testing fact sheets

http://www.cdc.gov/drugresistance/

**CDC** drug resistance


http://www.cdc.gov/drugresistance/community/

**CDC** antibiotic resistance

http://www.asmusa.org/division/c/index.htm

ASM Division C (clinical microbiology) website includes "Askit" feature

Thankyou



www.phppo.cdc.gov/dls/master/default.asp