

Abbreviated Identification of Bacteria and Yeast

NCCLS Guideline M 35 A

Mary K. York, Ph.D. ABMM

References Prove Value

- Doern, G., R. Vautour, M. Gaudet, and B. Levy. 1994. Clinical impact of rapid *in vitro* susceptibility testing and bacterial identification.
 - ◆ MIC 9.6 h from colonies vs. 25.9 h- showed less mortality, length of stay and orders of laboratory tests
- Barenfanger, J., C. Drake, and G. Kacich. 1999. Clinical and financial benefits of rapid bacterial identification and antimicrobial susceptibility testing.
 - ◆ evaluated evening vs. next day - 5 h difference - length of stay and cost was significantly less

Pros and Cons of Rapid Methods

→ Pros

- ◆ Less work than standard methods
- ◆ Results are out faster
- ◆ Less cost

→ Cons

- ◆ Requires technical expertise for accuracy
- ◆ Cannot be applied to all cases
- ◆ May disrupt workflow

Criteria Used by Committee for Rapid Tests

- Only for specific organisms
- Errors must not have a negative impact on patient care
- Accuracy must be $> 95\%$
- Emphasis on organisms that have unique reactions
- Results are not presumptive if all criteria met-
cpt 4 code issue

Factors to Keep in Mind

- Not for direct specimens
- All conditions must be met
- Keep isolate for future testing if needed

Technologist must....

- Begin with pure colony
- Recognize what it could be from typical colony morphology
- Perform rapid tests accurately and read them correctly
- Often do Gram stain

Supervisor must.....

- Validate the competency of the staff doing tests
- Check to see that all tests are done
- Be sure that procedures are written and QC is done at appropriate intervals

CLIA '88

- *All quality control activities must be documented.*
- *The laboratory must check positive and negative reactivity with control organisms*
- *Each new lot/shipment of reagents, commercial tests, or biochemical test media prior to being used on patient specimens.*

CLIA '88

→ *and*

(1) Each day of use for catalase, coagulase, beta-lactamase, and oxidase reagents and DNA probes;

(2) Each week of use for Gram stain, bacitracin, optochin, ONPG, X and V discs or strips; and

(3) Each month of use for antisera...

Does not address ID disks, rapid indole, Staph Latex reagents, etc.....

CLIA 2003

→ *All quality control activities must be documented*

♦ *The laboratory must check positive and negative reactivity with control organisms*

(1) *Each day of use for DNA probes and beta-lactamase (ex. cefinase);*

(2) *Each week of use for Gram stain and AFB stains; and*

(3) *Every 6 months of use for antisera...*

(4) *Each new lot/shipment of reagents, commercial tests, or biochemical test media prior to being used on patient specimens.*

Organisms covered by NCCLS

Gram Negative Bacilli

- *Escherichia coli*
- *Haemophilus influenza*
- *Moraxella catarrhalis*
- *Proteus mirabilis/penneri*
- *Proteus vulgaris*
- *Pseudomonas aeruginosa*

Gram Positive Cocci

- *Enterococcus* species
- *Staphylococcus aureus*
- *Streptococcus agalactiae*
- *Streptococcus pneumoniae*
- *Streptococcus pyogenes*

Yeast

- *Candida albicans*
- *Candida glabrata*
- *Cryptococcus neoformans*

Organisms covered by NCCLS

Anaerobic Gram Negative Bacilli

- *Bacteroides fragilis* group
- *Bacteroides urealyticus*
- *Bilophila wadsworthii*
- *Prevotella* species
- *Prevotella intermedia*
- *Porphyromonas* species
- *Fusobacterium nucleatum*

Anaerobic Gram Negative Cocci

- *Veillonella* species

Anaerobic Gram Positive Bacilli

- *Clostridium difficile*
- *Clostridium perfringens*
- *Clostridium septicum*
- *Clostridium sordellii*
- *Clostridium tetani*
- *Propionibacterium acnes*

Anaerobic Gram Positive Cocci


- *Peptostreptococcus* species




Start with indole and oxidase

Indole

Oxidase



Cinnamaldehyde reagent




Not as immediate

Disappears faster

Is the reagent required for anaerobic spot indole

Identification of *E. coli*




- Indole +
- Oxidase -
- Gram-negative rod
- Beta-hemolytic

= *E. coli*


Limitation: some *Proteus* and *Morganella* and all *Edwardsiella* are hemolytic. These species are lactose-negative.

What if it is not hemolytic?



Lactose +
and

PYR negative = *E. coli*




PYR reactions of GNR


| Genus | Indole | Lactose | PYR |
|---|--------|----------------|-----|
| <i>Citrobacter</i> | Some + | V | + |
| <i>E. coli</i> | + | V | - |
| <i>Klebsiella</i> | Some + | + | + |
| <i>Yersinia</i> | Some + | - ⁺ | + |
| <i>Enterobacter</i> | - | + | + |
| <i>Serratia</i> | - | + | + |
| <i>Morganella</i> | + | - | - |
| <i>Proteus</i> | Some + | - | - |
| <i>Providencia</i> | + | - | - |
| <i>Edwardsiella</i> (beta hemolytic) | + | - | - |
| <i>Shigella</i> | Some + | - | - |
| <i>Salmonella</i> | - | - ⁺ | - |



What if it is not lactose-positive?




CLSI Algorithm for E. coli: Indole +; Oxidase –



Algorithm for some labs...

Indole + and oxidase -



Limitations

- Beta hemolytic, lactose-negative should be limited to UTI and geographic areas that lack hemolytic *Morganella* and *Proteus*.
- Take colony from BAP that corresponds to MAC or EMB
- Do not use MUG for GI specimens, except to identify *E. coli* O157.

Accuracy of 1064 Indole + oxidase- strains = 99.7%

- 294 were hemolytic and *E. coli*.
- 628 were lactose positive and PYR negative and *E. coli*.
- 65 were MUG positive and *E. coli*.
- 13 were MUG negative and needed kit to identify as *E. coli*.
- 64 were not *E. coli* but 3 were called *E. coli*

Cost savings


- \$3100 in reagents - 987 kits omitted
- 70 hours of technologist time

York, M.K., E.J. Baron, M. Weinstein, R Thomson, and J.E. Clarridge. 2000. A Multi-Laboratory validation of rapid spot tests for identification of *Escherichia coli*. J Clin Microbiol 38: 3394-3398.

Spreading *Proteus*



MacConkey



BAP

Proteus identification=

Spreading colony plus....


- ◆ Indole-positive: → *Proteus vulgaris*
- ◆ Indole-negative: → Ampicillin-susceptible:
Proteus mirabilis
→ Ampicillin-resistant
 - Maltose-negative or Ornithine-positive: *P. mirabilis*
 - Maltose-positive or Ornithine-negative: *P. penneri*

Proteus species

| | Spreading on BAP | Indole | Ampicillin S | ODC | Maltose ferment |
|---------------------|---------------------|--------|--------------|-----|--------------------|
| <i>P. mirabilis</i> | 96 | 2 | 95 | 99 | 0 |
| <i>P. vulgaris</i> | 64 | 98 | 0 | 0 | 97 |
| <i>P. penneri</i> | 65-90 | 0 | 0 | 0 | 100 |

Rapid maltose

Positive



Rapid Urea



- Make suspension of organism
- Add disk
- Incubate 2 h
- Read urea
- Add HCl -not nec.
- Add FeCl_3


Other Rapid Urea Positive species

- *Brucella*
- *Helicobacter pylori*
- *Bordetella* species
- Some *Corynebacterium*



- Grows on BAP
- Coccobacillus
- Catalase -positive
- Oxidase - positive
- Urease - positive

- This is *Brucella* unless proven otherwise
- Work in a BSC
- Confirm with serology
- Reportable disease
- Bioterrorism



Brucella Diagnosis


- Symptoms are non-specific
- Onset is insidious
- Risks are eating raw dairy products or working in a microbiology laboratory
- Without diagnosis there is no appropriate treatment:

Doxycycline plus rifampin

Pseudomonas aeruginosa



- Oxidase +
- Indole -
- Metalic or mucoid
- Fruity odor




**CF patients:
Confirm with
colistin or
poly B disk**

Identification of *Haemophilus influenzae*

- Small gram-negative coccobacilli
- Growth of large colonies only on CHOC in 24 h or around staphylococci
- And.....



ALA test negative - 2 h 35°C read under UV light



1. *Francisella* grows on CHOC in 48 h but is small colony
2. Cannot differentiate *H. influenzae* from *H. haemolyticus*; the latter is hemolytic on horse or rabbit blood agar.

- You can either make the reagent or buy it
- Recipe in NCCLS M35-A and handout



Haemophilus species

| | SATELLITE V factor | ALA | LACT | UREA | IND | ODC | CAT |
|----------------------------|-----------------------|-----|------|------|-----|-----|-----|
| <i>H. influenzae</i> | + | - | - | V | V | V | + |
| <i>H. haemolyticus</i> | + | - | - | + | V | - | + |
| <i>H. parahaemolyticus</i> | + | + | - | + | - | - | + |
| <i>H. parainfluenzae</i> | + | + | - | V | V | V | V |
| <i>H. paraphrophilus</i> | + | + | + | - | - | - | - |
| <i>H. ducreyi</i> | - | - | - | - | - | - | - |

Horse Blood Agar



Case Study

- ◆ 36 y/o male with HIV
- ◆ Camped in Yosemite
- ◆ Non-healing, erythematous 3 mm cyst on neck
- ◆ FNA aspirate and biopsy
- ◆ Gram stain negative
- ◆ GNR grew on chocolate in 3 days
- ◆ Patient treated with ciprofloxacin & did well



Day 3 - GNR growing on Choc
catalase +; urease -; oxidase -
*Vitek NHI = 99% Actinobacillus
actinomycetemcomitans*


Day 4 - Satellite-negative
ALA neg; motility -; NH 2520 No ID

Day 5 - Tech sets up MIC - no growth day 7


Day 10 - KB set up

Day 12 - KB penicillin R

Day 40 - MDL reports ID-done by FA




Francisella tularensis - BSL III




Francisella tularensis

- ◆ If growth takes longer than 24 h, may be *Francisella* which will not satellite,
 - ALA-negative.
 - Grows on CHOC but not on BAP.
 - Oxidase-negative; catalase weak.
 - No kit will identify
 - Do cefinase-result is positive.
 - Send to health department if cefinase-positive
- ◆ Rabbit fever, tick, mosquito & fly bites
- ◆ Bioterrorism




Moraxella catarrhalis




Confirmatory Identification of *M. catarrhalis*

- Gram negative diplococci
- Grows on BAP
- Oxidase-positive
- Butyrate (or indoxyl acetate) positive in 5 min



Actually all *Moraxella* are butyrate + but the others are coccobacilli



- ◆ Inoculate plate and put Pen disk on it
- ◆ Perform Gram stain from around the disk
- ◆ Do on subs of positive blood cultures - polymyxin B and vancomycin too!



- May use indoxyl acetate
- Same method and color
- Can be used for identification of *Campylobacter jejuni*
(*C. lari* is negative)

Campylobacter jejuni

- ◆ Requires microaerobic environment
- ◆ Oxidase - positive
- ◆ Catalase - positive
- ◆ Curved rod
 - Hippurate - positive
 - If hippurate - negative; indoxyl acetate -positive and cefazolin R identifies
 - Nal Acid or cipro R & cefazolin R & indoxyl acetate -negative is *C. lari*



Rapid Hippurate


- Inoculate broth
- Incubate 2 h
- Add ninhydrin
- Observe blue color



In Summary for Gram negative rods

- If growing on MAC, do
 - ◆ Indole
 - ◆ Oxidase
 - ◆ Kit if not *E. coli*, *Proteus* or *Pseudomonas*
- If not growing on MAC, do
 - ◆ Catalase
 - ◆ Oxidase
 - ◆ Gram stain
 - ◆ Generally kits are not helpful here

Staphylococcus aureus



Latex agglutination

Identification of *S. aureus*


- Grows as opaque, white colony
- May be beta hemolytic
- Catalase-positive
- Gram-positive cocci in clusters
- Slide, tube or agglutination positive

Rapid tube coagulase is 4 h at 35°C; then must go to 25°C for 20 h more if negative

Limitations

- Rapid Tube Coag is 4 h at 35°C; then must go to 25°C for 20 h more if negative
- *S. lugdunensis* and *S. schleiferi* can be slide and Latex positive
- *S. saprophyticus* (and rare others) can be Latex positive
- MRSA can be Latex negative

Algorithm for Latex testing for Staphylococci



Hemolysis is defined as present in 18 h - not 48 h; not under colony

***Do PYR if tube coagulase negative to rule out *S. lugdunensis*⁵⁶**

Case study

- Positive blood culture from patient with endocarditis
- Slide or Latex coagulase positive
- Tube coagulase-negative
- PYR positive
- Ornithine positive

Staphylococcus lugdunensis

- Identification important to treatment
- CLSI now uses *S. aureus* methicillin breakpoints for this species





Use PBP2a on CoNS with MICs of less than 4 $\mu\text{g/ml}$



| | Hemolysis 18 h | Slide Coag | Latex agglu | Tube Coag | ODC | Poly B | PYR |
|---------------------------------|----------------|------------|-------------|-----------|-----|--------|-----|
| <i>S. aureus</i> | + | V | + | + | - | R | -? |
| <i>S. intermedius</i> (dogs) | V | V | V | + | - | S | + |
| <i>S. lugdunensis</i> | - | + | + | - | + | R | + |
| <i>S. schleiferi</i> | - | + | + | - | - | S | + |
| <i>S. saprophyticus</i> * | - | - | V | - | - | S | - |
| <i>S. epidermidis</i> | - | - | - | - | V | R | - |
| most others | - | - | - | - | - | S | V |


*Novobiocin resistant.

Streptococcus pneumoniae



30 min at 35°C 61

Quellung reaction



S. pneumoniae Identification

- Colonies small, transparent- may be mucoid
- Gram positive cocci - lancet shaped in pairs
- Catalase negative
- Bile soluble

Limitation: Not all are bile soluble but all bile soluble are *S. pneumoniae*


Enterococcus Identification

- Large (1mm) non hemolytic colony
- Gram-positive cocci in pairs and chains
- Catalase-negative
- PYR +



Case Study

- Positive blood culture post-partum
- PYR positive
- Laboratory called it *Enterococcus*
- Gram stain alerted to misidentification
- Called contaminant: patient discharged.



Other PYR-positive Cocci

Aerococcus - tetrads

Vagococcus - motile

Vancomycin resistant is always
Enterococcus if PYR +

Cannot separate *Enterococcus*
from *Lactococcus*




Phenotypic characteristics^a

| Genus | Gram | CAT | LAP | NaCl | 10°C | 45°C | Colony | Hemolysis on BAP |
|---|---------|-------|-----|------|------|------|------------------------|------------------|
| <i>Enterococcus</i> | ch | — | + | + | + | + | Large | α/γ/β |
| <i>Lactococcus</i> | ch | — | + | v | + | - | Large | α/γ |
| <i>Vagococcus</i> (motile) | ch | — | + | + | + | - | Large | α/γ |
| <i>Abiotrophia/Granulicatella</i> | ch | — | + | - | - | v | Satellite | α/γ |
| <i>Globicatella</i> | ch | — | - | + | - | - | Small | α |
| <i>Dolosicoccus</i> | ch | — | - | - | - | - | Small | α |
| <i>Aerococcus viridans</i> | cl/t | —,w | - | V | - | - | Large | |
| <i>Helcococcus kunzii</i> | cl/t | — | - | V | - | - | Small | γ |
| <i>Gemella</i> | cl/t/ch | — | v | - | - | - | Tiny, 48 h to grow | α/γ |
| <i>Facklamia</i> | cl/ch | — | + | + | - | - | Small | γ |
| <i>Alloiococcus</i> | cl/t | w,+ | + | + | - | - | Tiny, 72 h to grow | γ |
| <i>Ignavigranum</i> | cl/ch | — | + | + | - | - | Satellite (v) or small | γ |
| <i>Rothia mucilaginosa</i> ^b | cl | —,w,+ | + | - | | | Sticky | γ |
| <i>Dolosigranulum</i> | cl/t | — | + | + | - | - | Small | γ |


Identification of *S. pyogenes*

- Colonies: dry, peaked >0.5 mm in diameter
- Beta hemolytic
- Catalase-negative
- Gram-positive cocci in pairs and chains
- PYR-positive or positive particle agglutination

Enterococcus can be beta




Enterococcus look different



- *Enterococcus* colonies are larger with less defined β zone
- Do rapid esculin if there is a concern
- *Enterococci* are esculin-positive

Group A Strep

PYR - Positive Algorithm




Beta Streptococci

Group B




Group G




- Catalase-negative
- Gram-positive cocci in pairs and chains
- Characteristic colony

and CAMP test positive

CAMP +



Reverse
CAMP +




Rapid CAMP-20 minutes at 35°C

...or hippurate positive


**Limitation: Do not
do hippurate on non-
hemolytic colonies-
Other Streptococcus
can be positive**



Candida albicans



Candida albicans




Only use 100% fetal or newborn calf serum-2 h 35°C

Candida albicans = *Positive for both enzymes*

Add developer
at 5 min

Read under
UV light
Like MUG

Use antibiotic
free media




Candida albicans Limitations

- Always start with a wet mount to confirm presence of yeast
- Cannot separate from *C. dubliniensis* which does not grow at 42°C-Not usually necessary
- *C. tropicalis* can form fringe (not feet) or get projections in germ tube after 3 h.

Candida glabrata

- Colonies smaller on BAP than other yeast
- Yeast in wet mount are tiny with no hyphae
- Colonies are larger on EMB than BAP at 24 h or....


Represents 20% of yeast in urine



Growth in 24 h at 35°C


Candida glabrata

- Colonies smaller on BAP than other yeast
- Yeast in wet mount are tiny
- Positive (yellow) RAT test in 3 h at 42°C



Cryptococcus neoformans

- Large mucoid colonies
- No pigment
- Capsule by India Ink or no pseudohyphae
- Round cells
- And positive caffeic acid test



Cryptococcus neoformans



Overnight



- Rapid caffeic acid
- Use media without dextrose
- 30°C 4 h is best

L-DOPA disk works better



Rapid testing can make a
difference

Thank you